



**CLAREON® PANOPTIX®
AND CLAREON®
PANOPTIX® TORIC**



EXECUTIVE SUMMARY



PanOptix® Trifocal IOL provides distance, 60 cm intermediate, and 40 cm near focal points¹



With **PanOptix®**, 20/20 visual acuity is possible at distance, intermediate, and near²



91.6% of PanOptix® recipients never required spectacles at any distance³



Clareon® IOLs have exceptional rotational stability, reduced edge glare, low PCO^{4,5} and are glistening free^{*6,7}



PanOptix® features a 60 cm intermediate focal point,^{1,8} ideally suited to tasks such as working on a computer⁹



99% of PanOptix® recipients would have the same lens implanted again¹

*Defined as modified Miyata grade 0, <25mv/mm² over 3 years (n=138), and over 9 years (n=20), respectively. †Based on n=1291 PanOptix® recipients. MVs, microvacuoles; PCO, posterior capsule opacification.

1. Modi S et al. *Ophthalmology* 2021;128:197; 2. Kohnen T, et al. *Clin Ophthalmol* 2023;17:155-163. 3. Zhu D, et al. *Ophthalmol Ther*. 2023; Published online ahead of print. doi: 10.1007/s40123-023-00657-5. 4. Alcon Data on File, 2019; 5. Lehmann R et al. *Clin Ophthalmol* 2021;15:1647; 6. Stanojic N et al. *J Cataract Refract Surg* 2020;46:986; 7. Oshika T J et al. *Cataract Refract Surg* 2020;46:682; 8. Lwowski C et al. *J Cataract Refract Surg* 2021. doi:10.1097/j.jcrs.0000000000000780; 9. Charness N et al. *Proc Hum Factors Ergon Soc Annu Meet* 2008;52:1614. 10. Alcon Data on File, 2022.

A woman is sitting on a dark, rocky cliff edge, looking out over the ocean at sunset. The sun is low on the horizon, creating a warm, golden glow across the sky and reflecting on the water. The sky is filled with soft, wispy clouds. The woman is wearing a light-colored t-shirt and shorts. The overall mood is contemplative and serene.

PATIENT AND ECONOMIC BURDEN OF PRESBYOPIA


PRESBYOPIA IS AN AGE-RELATED IMPAIRMENT

of near vision characterized by a gradual decrease in accommodation of the eye

Clear vision



Presbyopia



Distance (∞)

Intermediate (60 cm)

Near (40 cm)




Presbyopia reduces the eye's ability to stretch the lens and focus on near objects, as a result of age-related changes to the natural elasticity of the lens




UNCORRECTED PRESBYOPIA


can reduce a patient's quality of life and impacts day to day activities




Feel **disabled** by reduced near vision*¹



Say presbyopia **impacts** their **ability to read**^{†2}




Say presbyopia **impacts** their **ability to use digital devices**^{†2}



Report a **decrease in quality-of-life** score*¹



Report a **low sense of accomplishment** due to impaired vision*¹




Require help in performing routine activities of daily living*¹

**CLAREON® PANOPTIX®
OPTICAL DESIGN**



TRADITIONAL TRIFOCALS


use sequential diffractive optics that require an 80 cm intermediate focal point for a 40 cm near focal point



It has been demonstrated that there is better comfort at 60 cm for intermediate tasks than at 80 cm²⁻⁵

A TRIFOCAL USING SEQUENTIAL DIFFRACTIVE OPTICS


can obtain a more comfortable intermediate focal point of 60 cm with a near focal point of 30 cm



A near focal point of 33 cm is not recognized by major organizations such as ANSI and AOA, which both define near vision as 40 cm^{2,3}

CREATING A FOURTH FOCAL POINT


at 120 cm enables a new 60 cm intermediate, while maintaining a 40cm near focal point



However, energy being directed to a fourth focal point could compromise distance visual performance



PROPRIETARY ENLIGHTEN® OPTICAL TECHNOLOGY

redirects light energy from the 120 cm focal point to distance¹



This results in the Clareon® PanOptix® Trifocal IOL design with distance, 60 cm intermediate, and 40 cm near focal points¹

PANOPTIX® TRIFOCAL IOL IS A SINGLE-PIECE, ASPHERIC, NON-APODIZED DIFFRACTIVE LENS



PanOptix® has a 6.0-mm biconvex optic, a central trifocal zone of 4.5 mm, an overall diameter of 13.0 mm, and 0° haptic angulation

PANOPTIX® USES THE PROPRIETARY ENLIGHTEN® TECHNOLOGY¹



PanOptix® transmits 88% of available light to the retina,* 50% of which is allocated to distance vision, 25% to intermediate vision, and 25% to near vision^{1,2}


*At the simulated 3.0-mm pupil size.¹

1. Modi S et al. *Ophthalmology* 2021;128:197; 2. Alcon Data on File, 2015.

CLINICAL DATA



A POOLED ANALYSIS OF MULTIPLE STUDIES DEMONSTRATED THAT PANOPTIX® PROVIDES A CONTINUOUS RANGE OF VISION FROM DISTANCE TO NEAR




20/20 visual acuity is now possible at distance, intermediate (66cm), and near (40cm)[†]

*Pooled analysis data is from six individual studies¹⁻⁶, detailed in the slide notes; †Snellen VA was converted from logMAR VA. A Snellen notation of 20/20-2 or better indicates a logMAR VA of 0.04 or better, which means ≥ 3 of the 5 ETDRS chart letters in the line were identified correctly. logMAR, logarithm of the minimum angle of resolution; VA, visual acuity.

1. Kohnen T, et al. *Clin Ophthalmol* 2023;17:155-163.

PANOPTIX® DELIVERS EXCELLENT VISION AT -3.0 D (33CM)




PanOptix® demonstrates 20/25 (0.1 logMAR) VA at -3.0 D (33 cm) across multiple studies

*Pooled analysis data is from six individual studies¹⁻⁶, detailed in the slide notes; †Snellen VA was converted from logMAR VA. A Snellen notation of 20/20-2 or better indicates a logMAR VA of 0.04 or better, which means ≥ 3 of the 5 ETDRS chart letters in the line were identified correctly. logMAR, logarithm of the minimum angle of resolution; VA, visual acuity.

1. Kohnen T, et al. *Clin Ophthalmol* 2023;17:155-163.

PANOPTIX® DEMONSTRATES SUPERIOR VA AT INTERMEDIATE DISTANCES COMPARED WITH AT LISA* AND FINEVISION*



PanOptix® maintains a CDVA of 20/25 (0.1 logMAR) or better between 66 cm (-1.50 D) and 50 cm (-2.00 D), which is significantly better compared with AT LISA* and FineVision*


*Trademarks are the property of their respective owners.

CDVA, corrected distance visual acuity; logMAR, logarithm of the minimum angle of resolution; VA, visual acuity.

Carneros-Llorente AM *et al.* J Cataract Refract Surg 2019;45:587.

PANOPTIX® HAS SIMILAR BINOCULAR PHOTOPIC CONTRAST SENSITIVITY COMPARED WITH A MONOFOCAL CONTROL

Mean ($\pm 95\%$ CI) photopic binocular contrast sensitivity at Month 6*¹





There was no clinically meaningful difference in photopic contrast sensitivity between the two groups across spatial frequencies, regardless of the presence of glare^{†1}

*Best-case analysis set; †Losses of >0.3 log units at two or more spatial frequencies are considered clinically meaningful.² CI, confidence interval; cpd, cycles per degree. 1. Modi S et al. *Ophthalmology* 2021;128:197; 2. International Organization for Standardization (ISO). 11979-7:2014. Ophthalmic implants--Intraocular lenses--Part 7: Clinical investigations of intraocular lenses for the correction of aphakia. Available at <https://www.iso.org/standard/55684.html> (accessed 11 Jan 2022).

PANOPTIX® HAS SIMILAR BINOCULAR MESOPIC CONTRAST SENSITIVITY COMPARED WITH A MONOFOCAL CONTROL

Mean ($\pm 95\%$ CI) mesopic binocular contrast sensitivity at Month 6*








There was no clinically meaningful difference in mesopic contrast sensitivity between the two groups across spatial frequencies, regardless of the presence of glare^{†1}

ASSESS VISUAL DISTURBANCES WITH PANOPTIX®¹

A questionnaire was developed and validated to assess visual disturbances with PanOptix®¹. Questionnaire for Visual Disturbance (QUVID) contains **picture-referenced items** for patients to report the **frequency**, **severity**, and **bothersomeness** of glare, halo, starburst, hazy vision, blurred vision, double vision, and negative dysphotopsia²

3.1 In the past 7 days, <u>how often</u> did you experience glare?				
0	1	2	3	4
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Never	Rarely	Sometimes	Most of the time	Always


3.2 Using the pictures below, rate how severe your <u>worst experience</u> was with glare in the past 7 days. These pictures may not look exactly like what you see, but are a guide to help you choose your answer.				
0	1	2	3	4
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
None	A little	Mild	Moderate	Severe
				

3.3 In the past 7 days, how much has your glare <u>bothered</u> you?				
0	1	2	3	4
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Not bothered at all	Bothered a little bit	Bothered somewhat	Bothered quite a bit	Bothered very much

Patients can rate the frequency, severity, or bothersomeness of a particular visual disturbance

SEVERE VISUAL DISTURBANCES WERE EXPERIENCED BY LESS THAN 1 IN 6 PANOPTIX[®] PATIENTS IN THE FDA REGISTRATION TRIAL


Percentage of PanOptix[®] patients with a “severe” visual disturbance at Month 6*



≤16% of PanOptix[®] patients experienced a “severe” visual disturbance

LESS THAN 4.8% OF PANOPTIX® PATIENTS WERE VERY BOTHERED BY ANY VISUAL DISTURBANCES IN THE FDA REGISTRATION TRIAL


Percentage of PanOptix® patients “bothered very much” by visual disturbances at Month 6*




More than 95% of PanOptix® patients were not “bothered very much” by any visual disturbance

SEVERITY AND BOTHER FROM HALO PROFILES ARE MORE FAVORABLE WITH PANOPTIX® COMPARED WITH AT LISA*

Percentage of patients reporting “not at all” for halo severity



Percentage of patients reporting “not at all” for halo bother



Compared with AT LISA* patients, 19.7% more PanOptix® patients did not experience any “severe” halos and 17.4% more PanOptix® patients were not bothered at all by any halos

*Trademarks are the property of their respective owners.

†Pearson’s chi-squared test; ‡Mantel-Haenszel test.

Asena BS. *J Cataract Refract Surg* 2019;45:587.

NO HAZE, BLUR, DOUBLE VISION, OR DARK AREAS

No PanOptix® patients reported severe or very bothersome hazy vision, blurred vision, double vision, or dark areas

0.0%

Hazy vision, blurred vision, double vision, or dark areas

Percentage of **PanOptix® patients** with a “severe” visual disturbance at Month 6*

0.0%

Hazy vision, blurred vision, double vision, or dark areas

Percentage of **PanOptix® patients** “bothered very much” by a visual disturbance at Month 6*

0% of patients implanted with PanOptix® in the US FDA Registration trial reported “severe” or “very bothersome” hazy vision, blurred vision, double vision, or dark areas

PANOPTIX® PATIENTS CONSISTENTLY DEMONSTRATE SPECTACLE INDEPENDENCE AT ALL DISTANCES¹

In a meta-analysis of 13 unique clinical studies, 91.6% of PanOptix® IOL patients achieved complete spectacle independence.¹



PanOptix®

IOL is the most studied trifocal IOL technology¹



13

studies



513

patients from
14 countries



91.6%

of patients achieved
spectacle independence
across all distances

91.6% of PanOptix® patients reported never requiring spectacles for any distance¹

ASSESS PATIENT EXPERIENCE AND SATISFACTION WITH PANOPTIX®¹

A questionnaire was developed and validated to assess patient experience and satisfaction with panoptix®¹. The Intraocular Lens Satisfaction (IOLSAT) questionnaire **assesses a patient's subjective quality of vision** with respect to need for spectacles, as well as expectation and satisfaction with their quality of vision²

21. Overall, in the past 7 days, how satisfied were you with your vision?				
0	1	2	3	4
<input type="checkbox"/> Very dissatisfied	<input type="checkbox"/> Dissatisfied	<input type="checkbox"/> Neither satisfied nor dissatisfied	<input type="checkbox"/> Satisfied	<input type="checkbox"/> Very Satisfied
22. Given your vision today, if you had to do it all over, would you have the same lenses implanted again?				
0		1		
<input type="checkbox"/> No		<input type="checkbox"/> Yes		
23. Given your vision today, would you recommend the lenses you had implanted to your family or friends?				
0		1		
<input type="checkbox"/> No		<input type="checkbox"/> Yes		

Patients can rate how satisfied they are with their vision


PATIENT SATISFACTION IS HIGH AMONG PANOPTIX® PATIENTS¹

Delight your patients with the most implanted trifocal IOL worldwide²



95.3%

were satisfied with their vision^{*†‡}



99%

said they would have the same lens implanted again^{*†‡}



98.4%

said they would recommend the lens to family or friends^{*†‡}

The majority of PanOptix® recipients are satisfied with their vision, would have the same lens implanted again, and would recommend the lens to their family or friends (n=127)

*Month 6 all-implanted analysis set; †Assessed using the validated IOLSAT questionnaire; ‡n=127.
FDA, US Food and Drug Administration; IOLSAT, intraocular lens satisfaction.

1. Modi S et al. *Ophthalmology* 2021;128:197.
2. Alcon Data on File, 2022.

DELIVER EXCEPTIONAL VISUAL OUTCOMES WITH CONFIDENCE¹

PanOptix® is the most implanted trifocal IOL worldwide.³

Provide your patients with²⁻⁴:

- Exceptional patient satisfaction
- High rates of spectacle independence.




1. Kohnen T, et al. *Clin Ophthalmol* 2023;17:155-163.

2. Clareon® PanOptix® Directions for Use.

3. Alcon Data on File, 2022.

4. Zhu D, et al. *Ophthalmol Ther.* 2023; Published online ahead of print. doi: 10.1007/s40123-023-00657-5


CLAREON[®] MATERIAL AND DESIGN




CLAREON® MATERIAL AND DESIGN

- Aspheric hydrophobic acrylic IOL¹
- Ultraviolet and blue light filtering hydrophobic acrylate/methacrylate copolymer¹
- Index of refraction 1.55 at 35°C¹
- STABLEFORCE™ modified-L haptics¹
- Low levels of surface haze²
- Glistening free^{*3,4}

* Defined as modified Miyata grade 0, <25mv/mm² over 3 years (n=138), and over 9 years (n=20), respectively.
1. Clareon® PanOptix®. Directions for use. Available at www.ifu.alcon.com (accessed 11 Jan 2022); 2. Werner L et al. *J Cataract Refract Surg* 2019;45:1490; 3. Stanojic N et al. *J Cataract Refract Surg* 2020;46:986; 4. Oshika T et al. *J Cataract Refract Surg* 2020;46:682;



CLAREON® IOLS HAVE AMONG THE LOWEST LEVELS OF HAZE COMPARED WITH COMPETITOR IOLS



The ultra-smooth optic of Clareon® delivered significantly lower ($p < 0.001$) levels of surface haze compared with several competitor IOLs


*Surface haze and subsurface nanoglistenings were compared *in vitro*, n=10 IOLs per group.
CI, confidence interval.
Werner L et al. J Cataract Refract Surg 2019;45:1490.

CLAREON® IOLS DELIVER EXCEPTIONAL ROTATIONAL STABILITY THAT OPTIMIZES ASTIGMATISM CORRECTION

The high rotational stability of the Clareon® PanOptix® IOL instills the confidence that comes from predictable performance.^{1,2}

95.3% of Clareon® PanOptix® Toric IOLs rotate $\leq 5^\circ$ between surgery and day 1 measures.¹

Studies demonstrated an absolute rotation of **1.8°** at day 1 post-op and **2°** at 6 months post-op (n=124).²




1. Alcon Data on File, 2019.

2. Clareon® Toric IOL Directions for Use.

AXIAL DISPLACEMENT AND SIMULATED DIOPTRIC POWER SHIFT


Clareon® and AcrySof® Had the Lowest Axial Displacement and Simulated Dioptric Power Shift Compared with Competitor IOLs

Axial displacement and simulated dioptric power shifts for five compression diameters†




Clareon® IOL demonstrated significantly less axial displacement compared with enVista*, TECNIS*, and Vivinex* IOLs ($p < 0.001$)

CLAREON® IS DESIGNED TO REDUCE EDGE-ASSOCIATED GLARE AND PCO (1/2)



Proprietary edge curvature minimizes glare and positive dysphotopsias




Precision edge design guards against PCO and associated Nd:YAG procedures²

*Trademarks are the property of their respective owners.

†At a 550-illumination angle in vitro, n=5 IOLs per group. Nd:YAG, neodymium-doped yttrium aluminum garnet; PCO, posterior capsule opacification. 1. R. Nuijts. The Ophthalmologist. Clareon® IOL: A New Monofocal Platform. Highlights from Alcon's Satellite Symposium, held on October 9, 2017, at the XXXV Congress of the ESCRS, Lisbon, Portugal. 2017. Available at https://theophthalmologist.com/fileadmin/top/issues/1217/images/1217-900_Alcon_SS-new.pdf (accessed 11 Jan 2022); 2. Lehmann R et al. *Clin Ophthalmol* 2021;15:1647.


CLAREON® IS DESIGNED TO REDUCE EDGE-ASSOCIATED GLARE AND PCO (2/2)

Proprietary edge curvature minimizes glare and positive dysphotopsias



Precision edge design guards against PCO and associated Nd:YAG procedures

Clareon® Monofocal has shown low rates of PCO in a prospective, multicenter (16 sites), US study with 1 year of follow up:



SUMMARY



EXECUTIVE SUMMARY




PanOptix® Trifocal IOL provides distance, 60 cm intermediate, and 40 cm near focal points¹



With **PanOptix®**, 20/20 visual acuity is possible at distance, intermediate, and near²



91.6% of PanOptix® recipients never required spectacles at any distance³



Clareon® IOLs have exceptional rotational stability, reduced edge glare, low PCO^{4,5} and are glistening free^{*6,7}



PanOptix® features a 60 cm intermediate focal point,^{1,8} ideally suited to tasks such as working on a computer⁹



99% of PanOptix® recipients would have the same lens implanted again^{†1}

*Defined as modified Miyata grade 0, <25mv/mm² over 3 years (n=138), and over 9 years (n=20), respectively. †Based on n=1291 PanOptix® recipients. MVs, microvacuoles; PCO, posterior capsule opacification.


1. Modi S et al. *Ophthalmology* 2021;128:197; 2. Kohnen T, et al. *Clin Ophthalmol* 2023;17:155-163. 3. Zhu D, et al. *Ophthalmol Ther*. 2023; Published online ahead of print. doi: 10.1007/s40123-023-00657-5. 4. Alcon Data on File, 2019; 5. Lehmann R et al. *Clin Ophthalmol* 2021;15:1647; 6. Stanojic N et al. *J Cataract Refract Surg* 2020;46:986; 7. Oshika T J et al. *Cataract Refract Surg* 2020;46:682; 8. Lwowski C et al. *J Cataract Refract Surg* 2021. doi:10.1097/j.jcrs.0000000000000780; 9. Charness N et al. *Proc Hum Factors Ergon Soc Annu Meet* 2008;52:1614. 10. Alcon Data on File, 2022.

BACK-UP SLIDES



THE CUMULATIVE RATE OF SSIS WAS LOW AMONG PANOPTIX® RECIPIENTS



And similar to that of a monofocal control



The single SSI that occurred for the first eye for PanOptix® was an explant of the IOL due to subjective complaints of dissatisfaction with the level of vision

CLAREON® IOLS DELIVER AMONG THE LOWEST LEVELS OF SSNGS COMPARED WITH COMPETITOR IOLS

Comparison of glistenings in Clareon® and competitor IOLs



Clareon® IOLs demonstrated significantly less light scatter due to SSNGs compared with Eternity W-60, TECNIS OptiBlue, and EnVista IOLs

*Trademarks are the property of their respective owners.

N=10 per group.

CI, confidence interval; SSNG, subsurface nanoglistening. Werner L et al. J Cataract Refract Surg 2019;45:1490.

Alcon

SEE BRILLIANTLY